

**Product Insert**

**100bp DNA Ladder**


**Cat. No.**

EA23100 50ug (100Lane)

Shipping: On Dry/Blue Ice

Store at -20°C

**100bp DNA Ladder**



**Concentration: 100ng/µl**

**Recommended Loading**

3-5 µl/Lane.

**Concentration**

Typical Bands 100 ng/5µl

Other Bands 40 ng/5µl

**Recommended Electrophoresis Condition**

5 µl/Lane, 8 cm 1.7 %Agarose Gel, 0.5xTBE, 7V/cm, 1h.

**Contents (bp)**

100, 200, 300, 400, **500**, 600, 700, 800, 900, 1,000, 1,500

**Storage**

1 year at room temperature, for long-time preservation, please storage at -20°C .

**Description**

100bp DNA ladder is ideal for determining the size of double-stranded DNA from 100 to 1500 base pairs. The ladder consists of 11 linear double-stranded fragments .The **500**bp fragment is present at increased intensity to allow easy identification. All fragments are precisely quantified and mixed during the production. For 5 ul loading, all fragments except 500bp are 40ng. The 500bp fragment is 100ng. This ladder is pre-mixed with loading dye and is ready to use.

**Recommendations**

Do not heat before loading.

For quantification, adjust the concentration of the sample to equalize it approximately with the amount of DNA in the nearest band of the ladder. Visualize DNA by staining with ethidium bromide or with SYBR® Green I.

**Technical Support**

If the troubleshooting guide does not solve the difficulty you are experiencing, please contact Technical Support with details of reaction setup, cycling conditions and relevant date.

Address : Woolim Lions Valley A-606B, 168, Gasan digital 1-ro, Geumcheon-gu, Seoul, Korea 08507

E.mail : [info@philekorea.co.kr](mailto:info@philekorea.co.kr)

Website : [www.philekorea.co.kr](http://www.philekorea.co.kr)

Tel. +82 2 2105 7020 +82 42 862 9636

Fax. +82 2 2105 7025 +82 42 862 9638