

eCube Gel & PCR Purification Kit

\* for Gel Purification

eCube Gel & PCR Purification Kit

\* for PCR Purification

\* Things to do before starting

- For 50preps, add 60ml ethanol (96-100%) to Wash Buffer when first open.  
For 200preps, add 180ml ethanol (96-100%) to Wash Buffer when first open.  
For 300preps, add 280ml ethanol (96-100%) to Wash Buffer when first open.
- Prepare a 55°C dry bath or water bath.
- All centrifuge steps are done at full speed (**14,000 rpm or 10,000 x g**) in a microcentrifuge.

\* Things to do before starting

- For 50preps, add 60ml ethanol (96-100%) to Wash Buffer when first open.  
For 200preps, add 180ml ethanol (96-100%) to Wash Buffer when first open.  
For 300preps, add 280ml ethanol (96-100%) to Wash Buffer when first open.
- All centrifuge steps are done at full speed (**14,000 rpm or 10,000 x g**) in a microcentrifuge.

 1. Excise the **agarose gel (up to 300mg)** containing relevant DNA fragments with a clean scalpel.

↓

 2. Add **500ul of XPDF Buffer** to the sample and mix by vortexing.  
(For >2% agarose gels, add 1ml of XPDF Buffer.)

3. Incubate at **55°C for 10 - 15 min** and vortex the tube every 3 min until the gel slice dissolved completely.

↓

 4. Transfer the sample mixture to **XPDF Column**.


↓

 5. Centrifuge for **30 sec** then discard the flow-through

↓

 6. Add **750ul of Wash Buffer (ethanol added)** to the XPDF Column.


↓

 7. Centrifuge for **30 sec** then discard the flow-through

↓

 8. Centrifuge for **3 min** to dry column


↓


 9. Add **40ul of Elution Buffer or ddH<sub>2</sub>O** (pH 7.0 – 8.5) to the membrane center of XPDF Column. Stand the column for 2 min.

↓

 10. Centrifuge for **2 min**


↓

 11. You can get pure DNA.


 1. Transfer **up to 100ul of PCR product** (excluding oil) and add **5 volumes of XPDF Buffer** to a 1.5 ml microcentrifuge tube (not provided) then mix well by vortexing.

↓

 2. Transfer the sample mixture to **XPDF Column**.

 3. Centrifuge for **30 sec** then discard the flow-through


↓

 4. Add **750ul of Wash Buffer (ethanol added)** to the XPDF Column.


↓

 5. Centrifuge for **30 sec** then discard the flow-through

↓

 6. Centrifuge for **3 min** to dry column

↓

 7. Add **40ul of Elution Buffer or ddH<sub>2</sub>O** (pH 7.0 – 8.5) to the membrane center of XPDF Column. Stand the column for 2 min.

↓

 8. Centrifuge for **2 min**

↓

 9. You can get pure DNA.